
# Introcan Safety<sup>®</sup>

## Introcan Safety<sup>®</sup>-W

---

---

**B | BRAUN**

LLDorder 3846 - grunewald 171708

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




## GB Instructions for use

### Materials used

Introcan Safety® is available in two materials:  
FEP, PP, ABS, chrome-nickel steel  
PUR, PP, ABS, chrome-nickel steel

### Indications

The Introcan Safety® Intravascular catheter is inserted into a patient's vascular system for short term use (less than 30 days) to sample blood, monitor blood pressure, or administer fluids and blood intravascularly.

The Introcan Safety® is also indicated for subcutaneous infusion therapies.

14-24 gauge catheters may be used with power injectors for which the maximum pressure setting is 300 psi/21 bar.

Use luer lock connection only.

See table below for recommended maximum flow rates for the Introcan Safety® IV Catheter when using a power injector. Testing performed at room temperature (22°C / 72°F). Warming contrast media per manufacturers' recommendations, can reduce the pressure necessary to achieve recommended flow rates.

		Contrast Media (Viscosity at 20°C) [mPa*s]	Flow Rate (ml/sec)
24 Gauge	all lengths	2,3	5
		27,5	2,5
22 Gauge	G22x25	2,3	8
		27,5	5
	G22x64	2,3	7
		27,5	3
20 Gauge	G20x25	2,3	14
		27,5	10
	G20x32	2,3	13
		27,5	8
G20x45	2,3	19	
	27,5	15	
18 Gauge	all lengths	2,3	19
		27,5	15

The above are flow rates tested in a laboratory setting and are the maximum flow rates our IV Catheters are able to achieve and do not constitute guarantees, warranties, or predictions regarding the outcome of your matter.

In any case, it is the responsibility of the user to accommodate the flow rates to each patient's conditions and / or therapy needed.

Facilitate the placement of Vascular Access Devices such as guidewires, indwelling central venous catheters, peripherally inserted central catheters, and midline catheters into the vascular system.

### Contraindications

Introcan Safety® I. V. Catheter should not be used in patients with known hypersensitivity to any of the materials employed.

### Risks

This I. V. catheter is designed to reduce the risk of accidental needlesticks; however, care must be taken to avoid needlesticks. Report needlesticks immediately and follow established protocols. Standard precautions must be adhered to in accordance with Centers for Disease Control and Prevention / Occupational Safety and Health Administration (CDC/ OSHA) standards for bloodborne pathogens, when starting or maintaining any I. V. catheter, to avoid the risk of exposure to contaminated blood. Depending on how long the catheter is left in situ, on the type and amount of infusions or injections administered, and on individual predisposition, thrombophlebitis may occur in the accessed vein. In arterial puncture, the artery may occlude in rare cases due to thrombotic or embolic complications, resulting in ischaemia.

### Duration of use

Change according to CDC Guidelines and / or Hospital or Institutional protocols. The puncture site should be checked at regular intervals. Introcan Safety® I. V. Catheter should be removed in the event of local or systemic signs of infection.

### Warning

Re-use of single-use devices creates a potential risk of patient or user. It may lead to contamination and/or impairment of functional capability. Contamination and/or limited functionality of the device may lead to injury, illness or death of the patient.

In the unlikely event that the safety mechanism does not engage, keep the needlepoint away from the body and fingers at all times and immediately dispose into an approved sharps container.

After withdrawal, do not reintroduce the steel needle into the catheter, as the latter may be cut off, leading to catheter embolism.

Do not bend the catheter/needle during insertion, advancement, or removal of the needle. Extreme care should be taken not to cut the catheter and possibly cause an embolism. In the case of an unsuccessful IV start, remove the needle first to activate safety mechanism, then remove catheter from patient and discard both. Do not use scissors or sharp instruments at or near the insertion site.

Always clearly mark arterial lines to avoid inadvertent injection. Verify adequate collateral circulation prior to arterial puncture.

Luer-slip connections should not be left unattended due to potential of disconnection. If blood is present, flush device according to institutional protocol. Do not attempt to override or defeat the safety shield.

Prior to use with power injectors, ensure that secure connection exists between the catheter and power injector.

Measures should be taken to avoid kinking or obstructing the catheter system during power injection to avoid product failure.

Patency of catheter must be established prior to use with power injector.

Use only if packaging is intact.

### Guidelines

- This device is intended for single use only and is provided sterile and non-pyrogenic.

- This device is sterile unless package has been opened or damaged.
- Components do not contain natural rubber latex, PVC, or DEHP.

### Application

#### NOTE:

Do not rotate catheter prior to insertion.

1. After disinfection of the puncture site and removal of the protective guard, puncture a suitable vein. If venipuncture was successful, blood will immediately be visible inside the transparent flashback chamber (see figure A).
2. Lower and advance the entire catheter and needle unit slightly to ensure catheter tip is in vessel (see figure B).
3. Using the push-off plate, advance catheter off the needle (approximately 1/8" or 3mm) and observe blood return between the needle and catheter to confirm catheter is in the vessel (see figure C). After confirmation, continue advancing catheter off the needle and into vessel.
4. Release the tourniquet.
5. Before removing the steel needle, compress the vein distal to catheter tip to prevent spillage of blood. At the same time stabilize the catheter hub to prevent catheter dislodgement during needle removal. Withdraw needle straight back with a controlled and continuous motion (minimize rotation of needle). Metal safety shield will automatically attach to needle tip as needle tip exits catheter hub (see figure D).
6. Dispose of needle immediately into sharps container.
7. Connect to infusion line and cover puncture site with a sterile dressing (see figure E). If applicable, remove vented flashplug for glucose test.
8. If used to facilitate placement of a Vascular Access Device follow manufacturer's instructions for use of the specific device being placed.

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B. Braun Melsungen AG