


hCG Cassette (25 mIU/mL)

**Alere™** hCG Cassette (25 mIU/mL)**Intended Use**

**Alere™ hCG Cassette (25 mIU/mL)** is intended for the qualitative detection of human chorionic gonadotropin (hCG) in urine as an aid in the diagnosis of pregnancy. For professional *in vitro* diagnostic use only.

**Introduction**

hCG is a glycoprotein hormone produced by the blastocyst.<sup>1,2</sup> hCG normally begins to be detected in the urine from 7 days after conception. The sudden rapid rise in concentration of hCG in urine following conception makes it an excellent marker for pregnancy.<sup>3,4</sup>

**Test Principle**

**Alere™ hCG Cassette (25 mIU/mL)** is a rapid chromatographic immunoassay for the qualitative detection of human chorionic gonadotropin in urine to aid in the early detection of pregnancy. The test uses two lines to indicate results. The test line utilizes a combination of antibodies including a monoclonal hCG antibody to selectively detect elevated levels of hCG. The control line is composed of goat polyclonal antibodies and colloidal gold particles. The assay is conducted by adding a urine sample to the sample well of the test device and observing the formation of colored lines. The sample migrates via capillary action along the membrane to react with the colored conjugate.

Positive samples react with the specific antibody-hCG-colored conjugate to form a colored line at the test line region of the membrane. Absence of this colored line suggests a negative result. To serve as a procedural control, a colored line will always appear in the control line region indicating that the correct volume of sample has been added and membrane wicking has occurred.

**Kit Contents And Storage****Materials Provided**

Each **Alere™ hCG Cassette (25 mIU/mL)** kit contains sufficient materials for 20 tests: The device contains anti-hCG particles and anti-hCG coated on the membrane.

- 20 Pouched devices: Each sealed foil pouch contains 1 device, 1 disposable pipette and 1 desiccant packet
- 1 Package insert

Store at 2-30°C. Do not use after the expiry date.

**Materials Required But Not Provided**

- Sample collection container
- Clock, timer or stopwatch

**Precautions**

1. Do not open the foil pouch until ready to test.
2. Do not use devices that have become wet, or if the foil pouch has been damaged.
3. Do not use kit beyond expiration date printed on the outside of the kit carton.
4. Do not reuse device.
5. Do not eat, drink or smoke in the area where the specimens or kits are handled.
6. Handle all specimens as if they contain infectious agents. Observe established precautions against microbiological hazards throughout testing.
7. It is recommended that disposable gloves should be worn while handling specimens.

8. All specimens and specimen contaminated materials should be disposed of in accordance with local biohazard waste disposal protocol.
9. To obtain accurate results, the package insert instructions must be followed.

### Sample Collection And Storage

A urine sample collected at any time of the day is suitable, but a first morning sample is recommended.<sup>5</sup> Urine samples must be collected in a clean and dry container. Samples may be stored in the refrigerator (2-8°C) for up to 48 hours or frozen below -20°C. Samples must be allowed to reach room temperature (15-30°C) prior to testing. Urine samples exhibiting visible precipitates should be centrifuged, filtered, or allowed to settle to obtain a clear sample for testing. **Alere™ hCG Cassette (25 mIU/mL)** is not validated for use with samples containing preservatives.

### Assay Procedure

**Allow the device and urine sample to equilibrate to room temperature (15-30°C) prior to testing.**

1. Remove the device from the sealed pouch and place on a clean and level surface.
2. Using the plastic pipette supplied, insert the tip into the sample and squeeze the top bulb fully (A). DO NOT squeeze the bottom bulb (B) when using the pipette. Release the top bulb to draw up the liquid. The sample (100 µL) will be drawn into the lower part of the pipette. Excess sample will be drawn into the bottom bulb.
3. Remove the pipette from the sample. Squeeze the top bulb fully to dispense the sample **carefully** onto the sample well (S) and start the timer. Use a new pipette for each test performed, even if using the same urine sample.
4. **Read the result at 3 minutes after dispensing the sample. Do not interpret results after the 3 minute read time.** It is important that the background is clear before the result is read.

EN

### Interpretation Of Results

(Please refer to the illustration)

- **POSITIVE: Two distinct colored lines appear.** One line should be in the control line region (C) and another line should be in the test line region (T). The color intensities of lines may vary. Therefore, any shade of color in the test line region should be considered positive.
- **NEGATIVE: One colored line appears in the control line region (C).** No apparent colored line appears in the test line region (T).
- **INVALID: Control line fails to appear.** Insufficient sample volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new device. If the problem persists, discontinue using the kit immediately and contact your local distributor.

### Limitations

1. Positive results from very early pregnancy may later prove negative due to natural termination of the pregnancy.<sup>6</sup> It is therefore recommended that weak positive results are re-tested 48-72 hours later with a first morning urine sample.
2. A negative result may be obtained if a urine sample is too dilute. If pregnancy is still suspected, it is recommended the patient should be retested 48-72 hours later with a first morning sample.
3. Concentrations of hCG are generally lower in ectopic pregnancy than expected normal values for a given gestational age. Abnormal pregnancy cannot be distinguished from normal pregnancy by hCG levels alone.<sup>2,7</sup>
4. hCG remains elevated for a time after pregnancy.<sup>8</sup> Pregnancy tests carried out less than 3 weeks after giving birth or 9 weeks after natural loss or termination may need further evaluation.

5. A number of conditions other than pregnancy can cause elevated levels of urinary hCG e.g. menopause, trophoblastic disease and certain non-trophoblastic neoplasms.<sup>9</sup>
6. Samples containing <25 mIU/mL hCG may test positive, but samples containing <5 mIU/mL hCG should be negative.
7. Drugs containing hCG may interfere with **Alere™ hCG Cassette (25 mIU/mL)**, and produce misleading results.
8. False positive and false negative pregnancy tests may be observed in patients with abnormal bladder or kidney function e.g. enterocystoplasties<sup>10</sup> and renal failure.
9. If the test result with **Alere™ hCG Cassette (25 mIU/mL)** is not consistent with clinical evidence, further evaluation may be required.
10. Concentrations of hCG greater than 500,000 mIU/mL may elicit a prozone effect.
11. **Alere™ hCG Cassette (25 mIU/mL)** is not validated for use with samples containing preservatives.

## Expected Values

Urine samples from healthy males and post-menopausal females generally contain <10 mIU/mL hCG.<sup>11</sup> Levels are generally <5 mIU/mL in pre-menopausal females. On the first day of the first missed period, the levels of maternal urinary hCG are normally 50-250 mIU/mL. During the first trimester hCG levels peak at up to 200,000 mIU/mL in a typical pregnancy.<sup>3, 12</sup>

## Performance Characteristics

### Accuracy

A multi-center clinical evaluation was conducted comparing the results obtained using **Alere™ hCG Cassette (25 mIU/mL)** and another commercially available urine membrane hCG test. The urine study included 159 samples, and both assays identified 88 negative and 71 positive results. The results demonstrated a >99% overall

accuracy of **Alere™ hCG Cassette (25 mIU/mL)** when compared to the other urine membrane hCG test.

Method		Other hCG Rapid Test		Total Results
Alere™ hCG Cassette (25 mIU/mL)	Results	Positive	Negative	
	Positive	71	0	71
	Negative	0	88	88
Total Results		71	88	159

Sensitivity: >99.9% (95%-100%)\*

Specificity: >99.9% (96%-100%)\*

Accuracy: >99.9% (98%-100%)\*

\* 95% Confidence Intervals

### Sensitivity and Specificity

**Alere™ hCG Cassette (25 mIU/mL)** detects hCG at a concentration of 25 mIU/mL or greater. The test has been standardized to the W.H.O 4<sup>th</sup> International hCG Standard. The addition of LH (300 mIU/mL), FSH (1,000 mIU/mL), and TSH (1,000 µIU/mL) to negative (0 mIU/mL hCG) and positive (25 mIU/mL hCG) samples showed no cross-reactivity.

### Prozone Effect

**Alere™ hCG Cassette (25 mIU/mL)** has been shown to produce positive results with samples containing up to and including 500,000 mIU/mL hCG, which is higher than the maximum level expected during a typical pregnancy.

### Interfering Substances

The following potentially interfering substances were added to hCG negative and positive samples.

Acetaminophen	20 mg/dL	Caffeine	20 mg/dL
Acetylsalicylic Acid	20 mg/dL	Gentisic Acid	20 mg/dL
Ascorbic Acid	20 mg/dL	Glucose	2 g/dL
Atropine	20 mg/dL	Hemoglobin	1 mg/dL
Bilirubin	2 mg/dL		

None of the substances at the concentration tested interfered in the assay.

### Alere™ Product Support

Contact one of the following Alere™ Product Support Care Centers or your local distributor if you have any questions regarding the use of your Alere™ product. You may also contact us at [www.alere.com](http://www.alere.com).

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1. Hsu M., Kolm, P., Leete, J., Dong, K., Muasher, S., Oehninger, S. (1998). Analysis of Implantation in Assisted Reproduction through the use of Serial Human Chorionic Gonadotropin Measurements. *J. Assist. Reprod. Genet.* 15 (8): 496-503.
2. Catt K., Dufau M. & Vaitukaitis, J. (1975). Appearance of hCG in pregnancy plasma following the initiation of implantation of the blastocyst. *J. Clin. Endocrinol. Metab.* 40: 537-540.
3. Chard T. (1992). Pregnancy tests: review. *Hum. Reprod.* 7 (5): 701-710.
4. Corson S., Batzer, F. Schlaff, S. (1981). A Comparison of serial quantitative serum and urine tests in early pregnancy. *J. Reprod. Med.* 26 (12): 611-614.
5. Kaplan L.A. and Pesce A.J. (1989). The C. V. Mosby Company.
6. Wilcox A., et al. (1988). Incidence of early loss of pregnancy. *N. Engl. J. Med.* 319 (4): 189-194.
7. Braunstein G., Karow W., Gentry W., Rasor J., Wade M. (1978). First trimester chorionic gonadotropin measurements as an aid in the diagnosis of early pregnancy disorders. *Am. J. Obst. Gynec.* 131 (1): 25-32
8. Steier J. & Myking, O. (1984). Human Chorionic Gonadotropin in Maternal Plasma after Induced Abortion, Spontaneous Abortion, and Removed Ectopic Pregnancy. *Obstet. Gynecol.* 64 (3): 391-394.
9. Braunstein G., Vaitukaitis, J., Carbone, P., Ross, G. (1973). Ectopic Production of Human Chorionic Gonadotrophin by Neoplasms. *Ann. Intern. Med.* 78: 39-45.
10. Nethercliffe J., Trewick, A., Samuella, C., Leaver, R., Woodhouse, C. (2001). False-positive pregnancy tests in patients with enterocystoplasties. *BJU International.* 87 (9) : 780-782.
11. Alfthan H., Haglund, C., Dabek, J., Stenman, U. (1992). Concentrations of Human Choriogonadotropin, its beta-subunit and the core fragment of the beta-subunit in serum and urine of men and nonpregnant women. *Clin. Chem.* 38 (10) : 1981-1987.
12. Lenton, E., Neal, L., Sulaiman, R. (1982). Plasma concentrations of human gonadotropin from the time of implantation until the second week of pregnancy. *Fertil. Steril.* 37 (6): 773-778.

**REF** CV506788C

Abon Biopharm (Hangzhou) Co., Ltd.  
#198 12<sup>th</sup> Street East,  
Hangzhou Economic & Technological Development Area,  
Hangzhou, 310018, P.R.China  
[www.abon.com.cn](http://www.abon.com.cn)

MT Promedt Consulting GmbH  
Altenhofstrasse 80  
66386 St. Ingbert/Germany

© 2017 Alere. All rights reserved.  
The Alere Logo and Alere are trademarks of the Alere group of companies.

1156126701

Pantone 032

Pantone 032 25%

Pantone 7544 50%

K100

Attention/注意:

1.

By approving the enclosed design draft, you (ABON's Customer) accept all responsibility for the accuracy of the design. If an error is detected following the printing or manufacturing of a material, you (ABON's Customer) are responsible for the cost of any inventory which is deemed unsuitable for sale. 作为ABON的客户，一旦你批准所附的设计稿后，即代表你已接受承担设计正确性的所有责任。如物料在随后的印刷和生产过程中发现有任何不适合销售的错误，你将负责承担任何库存的费用。

2.

Regulatory owner must ensure the compliance with all applicable regulations in the distribution territory. The applicable regulation means all federal, state and local laws, ordinances, rules, regulations, and mandatory ISO standards applicable to the design, development, manufacturer, control and marketing of the Product in the Territory. 法规所有人必须确保其销售区域适用法规的符合性。适用法规指的是所有联邦/州/地方性法律、法令、条例、规章以及强制性的ISO标准，以用于产品在某区域的设计、开发、生产、控制和上市。

Regulatory owner's RA/RR: Regulatory owner's RA/Authorized Regulatory Representative

<input type="checkbox"/> US	<input checked="" type="checkbox"/> OUS	<input type="checkbox"/> China			
Description 描述	Alere ABON CE FHC-102 booklet for Ireland	Part Number PN 号码	1156126701	Size 尺寸	110x165mm
Printing Contents 印刷内容	/	L Number L 号码	/	Size 尺寸	/
Designer 设计者	Amy	Design Date/Version 设计日期/版本	July 02, 2017/B	Mold Num. 模具号	/
Artwork Checked By 设计审核		Material/ Checked By 材质/审核	70g双胶		
Approved by Regulatory owner's RA/RR / Date 法规所有人RA/RR 确认/日期		Approved by US Customer US 客户确认/日期			
Approved by Affiliates QA / Date 关联公司QA确认/日期		Approved by Marketing / Date 市场部确认/日期			
Approved by ABON RA / Date ABON 注册部 确认/日期		Approved by PMT / Date 产品管理确认/日期			
Approved By ABON QA/Date ABON QA确认/日期		Effective Date 生效日期			