

Eco Sac

The cost effective Tissue Retrieval System

DEPLOY

- using 5mm instruments

POSITION


- with our easy orientation tabs and black stitching to simplify identification of the Sac's mouth

RETRIEVE

- via laparoscopic incision reducing hospital stay, cost and complication




Safe and secure
closure with strong
monofilament
drawstring












One-piece tear
resistant design with
unique PUR coating
for maximum strength
and durability




Easy orientation
using coloured
tab



Eco Sac Features:

-  Comprehensive solution covering laparoscopic procedures from node biopsy to enlarged spleens.
-  Minimally invasive reducing length of stay, financial strain and improved patient experience.
-  Our unique coated ripstop nylon is resistant to tearing if punctured.
-  Includes X-ray opaque marker.
-  Easy orientation using coloured tab.
-  Polyurethane coating impermeable to biological fluids.
-  Safe closure with strong monofilament drawstring.
-  Simplified deployment via the mouth or base.
-  Comprehensive tissue capture with large mouth opening.



New Product Code	Trocar (mm)	Diameter (mm)	Length (mm)	Volume (ml)	Number Per Box	GTIN-13 Code
EMP50ECO-R	5	38	90	50	10	5060358120008
EMP60ECO-10	10	45	100	80	10	5060358121098
EMP70ECO-10	10	57	100	100	10	5060358120015
EMP100ECO-10	10	70	144	200	10	5060358120046
EMP110ECO-R-10	5	55	135	100	10	5060358120299
EMP140ECO-10	12	115	185	550	10	5060358120077
EMP145ECO	12	102	235	750	5	5060358120114
EMP170ECO	12	170	180	1200	5	5060358120145
EMP230ECO	12	165	340	3100	5	5060358120206
EMP600ECO-5	15	320	410	6000	5	5060358120282

To talk to Fannin about our cost effective Eco Sac range, please call us on **+44 (0)1275 878 801** or email **espiner@fannin.eu**
www.espinermedical.com